


NEW ORDER of DRAW when COLLECTING MULTIPLE TUBES



The recommended order of draw when collecting multiple blood samples tubes.

(1) Sterile tubes for cultures (2) Light blue top for citrate tubes, (3) Red top serum tubes & gold SST gel-separator (4) green top plasma gel-separator with heparin and green heparin tubes, (5) lavender top for EDTA tubes, (6) other additive tubes (i.e. gray top in diagram) depends on the manufacturer.

***If only drawing Light Blue top Vacutainer, you can draw a serum discard tube first.**

04/2004